

RADIOPHARMACY DETECTION EQUIPMENT

OBJECTIVES:

1. DESCRIBE THE GENERAL MAKEUP OF A GAS-FILLED DETECTOR
2. DISCUSS HOW A GAS-FILLED DETECTOR OPERATES
3. DESCRIBE THE TWO MAJOR REGIONS WHERE THE GAS-FILLED DETECTORS FUNCTION
4. EXPLAIN THE USAGE OF THE TWO RADIOPHARMACY GAS-FILLED DETECTORS
5. DESCRIBE THE STRUCTURE OF A DOSE CALIBRATOR
6. DISCUSS HOW TO MEASURE A RADIONUCLIDE

10 CFR 35.50 STATES: "A LICENSEE SHALL POSSESS AND USE A DOSE CALIBRATOR TO MEASURE THE ACTIVITY OF DOSAGES OF PHOTON-EMITTING RADIONUCLIDES PRIOR TO ADMINISTRATION TO EACH PATIENT OR HUMAN RESEARCH SUBJECT."

GAS-FILLED DETECTOR

MAKEUP:

SEALED CHAMBER ENCLOSURES GAS

POSITIVE AND NEGATIVE ELECTRODES

HIGH VOLTAGE POWER SUPPLY

CURRENT METER MEASURES RADIATION


PRINCIPLE OF OPERATION:

RADIATION ENTERS ENCLOSED GAS

RADIATION CAUSES IONIZATION

IONS MOVE TOWARD RESPECTIVE POLES

IONS MEASURED ON CURRENT METER



GAS DETECTOR RESPONSE TO HIGH VOLTAGE

IONIZATION REGION:


USED TO MEASURE LARGER AMOUNTS OF
RADIATION

SATURATION VOLTAGE APPLIED

ALL PRIMARY IONS COLLECTED

ALPHAS, BETAS, AND GAMMAS SEPARATED

DIFFERENT ENERGY GAMMAS SEPARATED



GEIGER-MUELLER REGION:

USED TO MEASURE SMALL AMOUNTS OF
RADIOACTIVITY

THRESHOLD VOLTAGE APPLIED

PRIMARY AND SECONDARY IONIZATION

GAS VOLUME COMPLETELY IONIZED

CAN'T SEPARATE TYPES AND ENERGIES OF
RADIATIONS ABSORBED

DOSE CALIBRATOR

- CYLINDRICALLY SHAPED, SEALED CHAMBER WITH A WELL CONTAINING GAS
- HIGH VOLTAGE APPLIED TO ELECTRODES
- ENERGY SPECIFIC SETTINGS
- ACTIVITY READOUT IN μCi , mCi , Ci

HOW TO USE A DOSE CALIBRATOR: TURN ON MAIN POWER

PLACE SYRINGE OR VIAL HOLDER IN WELL

SELECT APPROPRIATE SETTING

ZERO THE DOSE CALIBRATOR

MEASURE SYRINGE OR VIAL

READ ACTIVITY